


Increasing the Linear Dynamic Range of Sector Field ICP-MS

Meike Hamester

Thermo Electron Bremen, Germany


Benefits of ICP-MS

- Elemental determinations of almost the whole periodic table
 - *Multi-elemental analysis in a single analysis run*
 - *Lower detection limits compared to Graphite Furnace Atomic Absorption (GFAA) or ICP-OES*
 - *High productivity with > 30 samples/hour*
- Limitations due to spectral interferences
 - *From the matrix or argon plasma itself*
 - *Limits the accuracy*




The screenshot shows a software window titled "HUUB.met - Method Editor Finnigan ELEMENT". The window contains a periodic table where several elements are highlighted with red boxes. The highlighted elements are: H, Li, Be, Na, Mg, K, Ca, Sc, Ti, V, Cr, Mn, Fe, Co, Ni, Cu, Zn, Ga, Ge, As, Se, Br, Kr, Rb, Sr, Y, Zr, Nb, Mo, Tc, Ru, Rh, Pd, Ag, Cd, In, Sn, Sb, Te, I, Xe, Cs, Ba, La, Hf, Ta, W, Re, Os, Ir, Pt, Au, Hg, Tl, Pb, Bi, Po, At, Rn, Ce, Pr, Nd, Pm, Sm, Eu, Gd, Tb, Dy, Ho, Er, Tm, Yb, Lu, Th, Pa, U, Np, Pu, Am, Cm, Bk, Cf, Es, Fm, Md, No, and Lr.

Who uses ICP-MS ?



Argon Plasma as Ion Source

- Both ICP-MS techniques (Quadrupole and Magnetic Sector) use the same ion source
- Samples are analyzed as solutions or solids
- The Sample is introduced into an argon plasma as a fine aerosol, via a nebulizer and spray chamber or by laser ablation of solid samples
- Within the plasma the solvent is evaporated and the sample species are decomposed into their constituent atoms and ionised
- Ionisation process is extremely efficient in the plasma, and contributes to the high sensitivity of ICP-MS



- 27 MHz Argon Plasma
- 6000 K

The drawback of the Argon plasma as an ion source is the formation of spectral interferences

What are spectral interferences ?

- Molecular species which occur at the same nominal mass as the analyte, they are
 - *Formed in the plasma and/or interface*
 - *Matrix dependent*
 - *Plasma condition dependent*
- They originate from the Argon plasma and the sample matrix

Isotope	Interference
^{31}P	$^{15}\text{N}^{16}\text{O}$ $^{14}\text{N}^{16}\text{OH}$
^{44}Ca	$^{28}\text{Si}^{16}\text{O}$ $^{12}\text{C}^{16}\text{O}_2$ $^{14}\text{N}_2^{16}\text{O}$
^{56}Fe	$^{40}\text{Ar}^{16}\text{O}$ $^{40}\text{Ca}^{16}\text{O}$ $^{28}\text{Si}_2$
^{60}Ni	$^{44}\text{Ca}^{16}\text{O}$ $^{23}\text{Na}^{37}\text{Cl}$ $^{28}\text{Si}^{32}\text{S}$

Examples for spectral overlaps

The Biggest Problem in ICP-MS: Spectral Interferences

Instrumental Solutions


- *High Resolution (SF ICP-MS)*
- *Cold Plasma*
- *Collision Cell/Dynamic Reaction Cell (Q ICP-MS)*
- *Mathematical Calculation*

Sample Preparation (Time-consuming, error-prone)

- *Pre-concentration/matrix evaporation (off-line)*
- *Membrane desolvation (on-line)*
- *ETV*


Building Blocks of an ICP-MS

Inductively Coupled Plasma Mass Spectrometer



Thermo Electron Factory in Bremen, Germany

- Since 1948 producing mass spectrometers
- Organic (MAT 95), IRMS (Triton, Delta, MAT 253/271/281) LTQ-FTMS and since 1994 ICP-MS as single (ELEMENT) and multicollector (Neptune)
- 175 colleagues working on magnetic sector field technology



New Factory



Worldwide 350 ELEMENT



Key characteristics


- Interference free analysis
- Matrix independent
 - *physical separation of interferences*
 - *multielement settings*
- High Signal to Noise Ratios
 - *< 0.2 cps for all 3 resolutions*
 - *across the mass range*
- Sensitivity:
 - *LR: 3,000,000 cps per from ~ ⁴⁵Sc upwards*
 - *MR: 3,00,000 cps per ppb from ~ ⁴⁵Sc upwards*
 - *HR: 50,000 cps per ppb from ~ ⁴⁵Sc upwards*

Under same set of conditions


SF ICP-MS for Demanding Applications

- Simultaneous ultra-trace to matrix analysis
- Ultra-trace determinations in complex matrices
- Analysis of radionuclides
- Isotope ratio analysis
- Analysis of fast transient signals
 - laser ablation, HPLC
- High throughput analysis


Principle components of a SF ICP-MS



Dispersion and Focusing – Magnetic Sector



Electrostatic Sector – Dispersion and Angular Focusing






- Francis W. Aston
 - *Tandem Electric field – Magnet field*
 - *Resolution ~ 130 (1919)*
 - *Stable Isotopes*
- Arthur J. Dempster (1918)
 - *180° magnetic sector*
- Josef Mattauch / Richard Herzog (1934)
 - *Double focusing mass spectrometer*
- Alfred Nier / Edgar G. Johnson (1953)
 - *Double focusing mass spectrometer with corrected second order aberrations*


Physical, matrix independent separation of interferences

➤ Definite relation between analyte and interference

$R = m/\Delta M$




Isotope	Interference	Resolution
³¹ P	¹⁵ N ¹⁶ O	1457
	¹⁴ N ¹⁶ OH	967
⁴⁴ Ca	²⁸ Si ¹⁶ O	2687
	¹² C ¹⁶ O ₂	1281
	¹⁴ N ₂ ¹⁶ O	965
⁵⁶ Fe	⁴⁰ Ar ¹⁶ O	2503
	⁴⁰ Ca ¹⁶ O	2480
	²⁸ Si ₂	2957
⁶⁰ Ni	⁴⁴ Ca ¹⁶ O	3058
	²³ Na ³⁷ Cl	2409
	²⁸ Si ³² S	3292
	²⁵ Mg ³⁵ Cl	3467




R = > 300
> 4000
> 10000

^{60}Ni in Urine



Key: matrix independent analysis



Why Use High Resolution for ICP-MS ?

Example: analysis of chromium in blood

Chromium has four naturally occurring isotopes

Only ^{52}Cr and ^{53}Cr are available:

Element: Cr Ord. Number 24

Selected	Isotope:	Mass:	Abundance:
	Cr50	49.9455	4.3450
x	Cr52	51.9400	83.7890
x	Cr53	52.9401	9.5010
	Cr54	53.9383	2.3650

Select Isotope for Measurement

Q-ICP-MS presence of Cr in sample


53Cr to confirm

50Ti & 50V

54Fe



Chromium in Blood

Serum/Urine/Seawater/Soil/Plants.....



Resolution vs. Sensitivity


- The change in mass resolution is achieved by changing the width of the entrance and exit slits of the mass spectrometer:
 - The wider the slit the higher the sensitivity
 - Fixed ratio between resolutions \Rightarrow independent of mass and matrix



Signal to Noise Ratio Sector Field ICP-MS


100 ng/L Li, Y, In, Tl

(<0.2 cps in all resolutions)




Calibration Curve

- IAEA 421 (water), dilution 1+1, Ra-226, Standard Addition




Accuracy and precision CASS3 (0.3% saltmatrix)

“Dilute and shot”




Dual Mode Detection System




Mass independent cross calibration
= no cross calibration

Dual Mode Detection System





Finnigan ELEMENT 2 Dual Mode Detection System




Peak Top Hopping in Both Mode

- Transient signals
 - Laser ablation
 - Chromatography coupling
- Isotope ratio analysis

One channel per isotope





Determination of large isotope ratios




Peak Shape

- The change in mass resolution is achieved by changing the width of the entrance and exit slits of the mass spectrometer:
 - The wider the slit the higher the sensitivity
 - Fixed ratio between resolutions \Rightarrow independent of mass and matrix



Isotope Ratio Analysis

$^{233}\text{U}/^{238}\text{U}$ Ratio IRMM 072/8




Finnigan ELEMENT2: NIST 981 Isotope Ratio Analyses

- NIST981 (Natural Pb)
 - 1ppb Pb (1.6×10^6 cps)
 - 3 Measurements (each consisting of 10*1min analyses)
 - All Pb ratios measured
 - ^{204}Hg on ^{204}Pb corrected for on-line
 - No outlier rejection
 - No correction for mass bias

	Pb206/Pb204	Pb207/Pb204	Pb207/Pb206	Pb208/Pb206	Pb208/Pb204
Run 1	16.953611	15.541761	0.916724	2.163897	36.685854
Run 2	16.952012	15.539417	0.916673	2.164276	36.688759
Run 3	16.948784	15.539274	0.916837	2.163845	36.674535
Average	16.951469	15.540150	0.916745	2.164006	36.683049
%RSD	0.015	0.009	0.009	0.011	0.020
Reference NIST	16.937096	15.491345	0.914640	2.168100	36.721317
Accuracy %	-0.08	-0.32	-0.23	0.19	0.10

Extended Dynamic Range in the Finnigan ELEMENT XR




> ppq

> ppt


> ppb

> ppm


Finnigan ELEMENT XR Detection System




Automatic switch to Faraday Detector



'Triple' Detector Mode: ^{36}Ar (LR, log scale)




Finnigan ELEMENT XR : Triple Mode Detection system



Matrix and ultra-traces
in one analysis

- 0.2 cps to 1×10^{12} cps
 - complete analysis from 0.1 ppt to 0.1% in LR (solution)
- Minimum integration time:
 - counting: 0.1 ms
 - analog: 1 ms
 - Faraday: 1 ms
- *No decay time with Faraday detection system*
 - Due to integration circuit
- Automatic switching between detection modes
 - *no preset*
 - *< 1 ms to Faraday detection*
- Automatic cross calibration

Finnigan ELEMENT XR Cross Calibration



- Mass independent cross calibration
 - e.g. Argon
- Fast
 - 3 s for ACF and FCF determination
- Reliable
 - automatic
 - user independent

Finnigan ELEMENT2 Analysis Time

clinical - Method Editor Finnigan ELEMENT

File Edit Setup View Resolution Interference Correction Plot Help

H																			He
Li	Be									B	C	N	O	F	Ne				
Na	Mg									Al	Si	P	S	Cl	Ar				
K	Ca	Sc	Ti	V	Cr	Mn	Fe	Co	Ni	Cu	Zn	Ga	Ge	As	Se	Br	Kr		
Rb	Sr	Y	Zr	Nb	Mo	Tc	Ru	Rh	Pd	Ag	Cd	In	Sn	Sb	Te	I	Xe		
Cs	Ba	La	Hf	Ta	W	Re	Os	Ir	Pt	Au	Hg	Tl	Pb	Bi	Po	At	Rn		
Fr	Ra	Ac																	
			Ce	Pr	Nd	Pm	Sm	Eu	Gd	Tb	Dy	Ho	Er	Tm	Yb	Lu			
			Th	Pa	U	Np	Pu	Am	Cm	Bk	Cf	Es	Fm	Md	No	Lr			

LOW

18

clinical - Method Editor Finnigan ELEMENT

File Edit Setup View Resolution Interference Correction Plot Help

H																			He
Li	Be																		
Na	Mg																		
K	Ca	Sc	Ti	V	Cr	Mn	Fe	Co	Ni	Cu	Zn	Ga	Ge	As	Se	Br	Kr		
Rb	Sr	Y	Zr	Nb	Mo	Tc	Ru	Rh	Pd	Ag	Cd	In	Sn	Sb	Te	I	Xe		
Cs	Ba	La	Hf	Ta	W	Re	Os	Ir	Pt	Au	Hg	Tl	Pb	Bi	Po	At	Rn		
Fr	Ra	Ac																	
			Ce	Pr	Nd	Pm	Sm	Eu	Gd	Tb	Dy	Ho	Er	Tm	Yb	Lu			
			Th	Pa	U	Np	Pu	Am	Cm	Bk	Cf	Es	Fm	Md	No	Lr			

Medium

10

Setup Global Parameter

Runs and Passes | Acquisition Default Settings | Evaluation Default Settings

Mode: 1 2 3

Res. Switching Delay [s]: 1

Acquisition: Runs 3, Passes 3, Res. 3

Evaluation: Runs 3, Passes 3, Res. 3

Estimated Time

Time per Pass [min:sec:ms]	Time per Res. [h:min:sec]
Low Res.: 00 : 02 : 499	Low Res.: 00 : 00 : 22
Med. Res.: 00 : 02 : 090	Med. Res.: 00 : 00 : 19
High Res.: 00 : 01 : 741	High Res.: 00 : 00 : 16
Total [h:min:sec]: 00 : 00 : 59	

OK

clinical - Method Editor Finnigan ELEMENT

File Edit Setup View Resolution Interference Correction Plot Help



H																			He
Li	Be																		
Na	Mg																		
K	Ca	Sc	Ti	V	Cr	Mn	Fe	Co	Ni	Cu	Zn	Ga	Ge	As	Se	Br	Kr		
Rb	Sr	Y	Zr	Nb	Mo	Tc	Ru	Rh	Pd	Ag	Cd	In	Sn	Sb	Te	I	Xe		
Cs	Ba	La	Hf	Ta	W	Re	Os	Ir	Pt	Au	Hg	Tl	Pb	Bi	Po	At	Rn		
Fr	Ra	Ac																	
			Ce	Pr	Nd	Pm	Sm	Eu	Gd	Tb	Dy	Ho	Er	Tm	Yb	Lu			
			Th	Pa	U	Np	Pu	Am	Cm	Bk	Cf	Es	Fm	Md	No	Lr			

High

2

Element XR Detection System – Filter Lens



Filter Lens



Filter Lens: ICP at $^{40}\text{Ar}^{40}\text{Ar}$

-20 V

7 V




Applications of ELEMENT XR

- **Geological:**

- *Determination of majors, traces and ultra-traces in survey analyses, replacing complimentary analysis techniques (e.g. AAS or ICP-AES).*
- *Use of the matrix element as internal standard in laser ablation analysis:*
 - Na in fluid inclusions.
 - Al in melt inclusions.
 - Ca in bone / corals / otoliths etc.
 - C in diamond analyses
 - Large isotope ratios
- *Concentration determination in minerals by laser ablation.*
- *Elemental ratios by laser ablation (e.g. Ca / Sr etc).*

Finnigan ELEMENT XR with Laser Ablation



- **Finnigan ELEMENT XR**
 - *All of the advantages of the ELEMENT2:*
 - High sensitivity
 - High resolution for reliable interference removal
 - Dark noise independent of resolution
 - Fast scanning
 - *Faraday Detector*
 - Increased linear dynamic range ($>10^{12}$)
 - No decay time
 - 1 ms integration time
 - Automatic switching
 - Automatic cross calibration